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This product is intended for use in custom medical devices manufactured by licensed custom manufacturing orthopaedic technicians under their sole responsibility. Orliman does not guarantee the correct operation of the product if the instructions provided are not followed.

FITTING INSTRUCTIONS FOR THE ORTHOPEDIC TECHNICIAN

When fitting the product, you must adhere to the following instructions:

- 1-Take the insoles and retrocapital metatarsal bars out of the package.
- 2-Depending on the condition and the needs of the user, it may or may not be necessary to place the retrocapital metatarsal bars.
- 3-The retrocapital bars, due to the material they are made of, allow us to work them (cut, trim, polish, etc.). In the case of needing metatarsal domes instead of retrocapital bars, you could convert the latter into metatarsal domes, just as if you wanted to convert them into heel pads for more cushioning in the heel.
- 4-If you need to place the metatarsal bar/dome, mark on the template where you are going to place it and add glue in both areas.
- 5-Let the glue dry for a couple of minutes.
- 6-Then glue those metatarsal bars/domes to the insole.
- 7-Add glue to the rest of the insole and the upper lining.
- 8-Let the glue dry for a couple of minutes.
- 9-Lower the upper lining little by little from the front part of the insole towards the back part, without stretching the lining.
- 10-Let the glue act for a couple of minutes so that both parts adhere well.
- 11-Trim off the excess upper lining, leaving the upper lining about 3 mm wider than the base of the resin.
- 12-Finally, we will finish off all the edges with scissors or a polishing machine, adapting to the last of the shoe where the user will wear the insoles.
- 13-If a small modification in the shape of the insole is necessary, it can be done by heating the area of the insole with a heat gun at a temperature of 70-90°, through the bottom part of the insole, at a distance of about 5 cm without stopping moving the gun in order to avoid unnecessary overheating of the material and the consequent loss of its properties/characteristics. Do this until you feel that the area becomes conformable.

IMPORTANT: These data may vary depending on the brand, gun model and room temperature. Perform checks with the free hand to ensure that the insole acquires the desired plasticity.

WARNING: If you make these small modifications using a heat gun, it is important not to place it on the user until it has cooled down and has the desired shape.

WARNING: Avoid overheating the material, keeping the heat gun moving and at a distance of about 5 cm. Always heat the bottom part of the insole, to avoid heating the offloading areas and the upper lining and losing part of its properties.

WARNING: Avoid using the insoles without socks.

***The orthopaedic technician will ensure the correct positioning/adaptation of the product.**

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